

# UCSD-IU ALS Metabolomics Study

IRB# 2009735565

Phlebotomy and Sample Processing Orders

January 7, 2021

## Phlebotomy

1. 4 ml of venous blood drawn into a lithium-heparin vacutainer tube (green-top tube; without gel separator [use of a gel separator tube is not acceptable]; 13 x 75 mm; BD #367884; not sodium-heparin)
  - a. Gently invert x 10 to mix
  - b. Keep at room temperature until plasma separation. Do not refrigerate. Refrigeration can activate platelets.
  - c. Separate plasma and whole cells by centrifugation within 1 hour of whole blood collection
    - i. Plasma is for metabolomics
    - ii. Whole cells are for epigenetics
2. 2.5 ml of venous blood drawn into a PAXGene RNA tube kept in an upright position for at least 10 seconds for full collection during venipuncture (16 x 100 mm with 6.9 ml of preservative; BD #762165). See attached "PAXGene tube blood collection guide"
  - a. Gently invert x 10 to mix
  - b. Be sure the PAXGene tube is labeled with: patient name, date, time, and "Naviaux ALS Study"
  - c. The PAXGene tubes are for microRNA analysis

## Sample Processing

1. Green-top tubes for plasma and whole cells
  - a. Centrifuge for 10 minutes at 900 g (typically about 1200-1500 rpm for most clinical centrifuges). Do not centrifuge for longer than 10 minutes.
  - b. Using a sterile transfer pipet or Pipetman, transfer and divide the plasma into two, (not one) 1.8 ml Nunc cryotubes (Thermo, #375418) labeled "Heparin-plasma".
  - c. Keep the remaining red and white cell fraction, mix and transfer to cryotubes labeled "Heparin-cells"
  - d. Label the tubes with: patient name, date, time, and "Naviaux ALS Study"
2. PAXGene RNA tubes
  - a. No further processing is required

## Sample Storage

1. Labeled cryotubes (separate tubes for plasma and whole cells)
  - a. Heparin-plasma: store at -70° to -80°C until the total set of 120 patient samples (30 subjects x 4 samples/subject x 2 tubes/sample time = 240 cryotubes) has been collected. -20°C storage is not acceptable.
  - b. Heparin-cells: same as above

2. PAXGene RNA tubes
  - a. Store at -70° to -80°C until the total set of 120 PAXGene tubes (30 subjects x 4 samples/subject) has been collected

### Sample Shipment

1. When the completed set of samples is ready for shipment, please email Dr. Naviaux at [naviaux@health.ucsd.edu](mailto:naviaux@health.ucsd.edu) before shipping to alert Naviaux Lab personnel at UCSD and to coordinate shipping.
  - a. Plan to ship on a Monday or Tuesday to ensure the samples DO NOT arrive on a weekend.
    - i. **WARNING:** Weekend delivery can result in samples being stored in a non-refrigerated FEDEX warehouse for several days and can result in a complete loss of all samples.
2. After receiving email confirmation from the Naviaux Lab, ship the samples on at least 10 pounds of dry ice in a Styrofoam box with inner dimensions of at least 6" x 8" x 9" by FedEx overnight to:

Dr. Robert K. Naviaux, MD, PhD  
UCSD School of Medicine  
214 Dickinson St., Room C107  
San Diego, CA 92103  
Cell: 619-993-2904

FedEx account #1034-4986-3

Be sure to include "Room C107" in the address. Failure to do this will result in failed deliveries.